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essenbioscience.com/phagocytosis


Quick guide

1 CULTURE TARGET CELLS & INDUCE APOPTOSIS (OPTIONAL)


Target cells (e.g. Jurkat)을 cell culture 하고 cytotoxic agent (e.g. camptothecin)를 treat(Optional)하여 apoptosis 유도.

2 LABEL WITH INCUCYTE™ PHRODO® RED CELL LABELING KIT

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FOR RESEARCH USE ONLY. NOT FOR THERAPEUTIC OR DIAGNOSTIC USE.

IncuCyte™ pHrodo® Red Cell Labeling Kit (Cat No. 4649) Components	Component	Amount	Storage	Stability
IncuCyte™ pHrodo® Red Cell Labeling Dye	A	1 vial	-20 °C	6 months
Dimethyl sulfoxide (DMSO), anhydrous	B	150 µL	18 - 25 °C	6 months
IncuCyte™ pHrodo® Red Cell Wash Buffer	C	100 ml	18 - 25 °C	6 months
IncuCyte™ pHrodo® Red Cell Labeling Buffer	D	100 ml	18 - 25 °C	6 months

Product label licence

RESEARCH Field This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not sell or otherwise transfer this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or

screening services, or information in return for compensation on a per-test basis; (c) manufacturing or quality assurance or quality control, or (d) resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.