

IncuCyte™ Jurkat NuLight™ Red Cells

Catalog number: 4613

Immortalized Human T lymphocyte cells stably expressing nuclear restricted red fluorescent protein.

Presentation, storage and stability

Jurkat NuLight™ Red Cells are supplied as 1 mL cryopreserved vials (1.5×10^6 cells/mL in 90% FBS and 10% DMSO) containing a stable population of human T lymphocyte cells expressing the NuLight™ Red fluorescent protein, restricted to the nucleus. NuLight™ cells should be stored in liquid nitrogen (vapour phase). Do not store at -80°C . When stored as recommended, the NuLight™ cells will be viable for at least 1 year from the date of receipt.

Background and intended use

Parental Jurkat (clone E6.1) cells (ATCC, Cat# TIB-152) were transduced with the IncuCyte™ NuLight™ Red Lentivirus Reagent (Cat# 4476; EF1 α , puromycin). 48 hours after transduction, the cell population was grown in complete growth media containing 1 $\mu\text{g}/\text{ml}$ puromycin, for 3 to 5 days, to select for cells expressing NuLight™ Red. Following selection the Jurkat NuLight™ Red Cells were validated by comparing morphology and chemotactic migration profiles to the parental cell line. The Jurkat NuLight™ Red Cells are certified mycoplasma free and are STR authenticated by ATCC.

Recommended use


Jurkat NuLight™ Cells should be thawed quickly in a 37°C water bath. Jurkat NuLight™ cells should be maintained at a density of between $0.1\text{-}1 \times 10^6$ cells/ml in growth medium containing $0.5 \mu\text{g}/\text{mL}$ puromycin. These cells are fully validated for use with the IncuCyte ZOOM® live-cell imaging and analysis system and are ideal for use in real-time cell counting studies, co-cultures and chemotaxis assays.

Please see the relevant protocol published on our website:
essenbioscience.com/nuilight

Safety data sheet (SDS) information

The SDS can be found on our website:
essenbioscience.com/nuilight

Parent Jurkat



Jurkat NuLight™ Red Cells







Figure 1. Representative images of Jurkat parent and Jurkat NuLight™ Red cells (fluorescence/phase blend and fluorescence only) and growth characteristics. Cells were seeded in PLO (0.01%) coated flat bottom plates and images were acquired using an IncuCyte ZOOM® system with 10x objective every 4 h. Growth curves were constructed from confluence data for both cell lines and red object count for Jurkat NuLight™ Red cells.




Recommended Media and Components

- RPMI (Cat# 21875 (UK) 11875 (US) Gibco/Life Technologies)
- 10% FBS (Cat# SH30071 Thermo Hyclone)
- 0.5 µg/ml Puromycin (Cat# A11138-03 Gibco/Life Technologies)

Growth Media

- RPMI + 10% FBS + 0.5µg/ml Puromycin

A. Jurkat Parent



B. Jurkat NuLight™ Red




Figure 2. Expression of NuLight™ Red has no impact on the chemotactic profiles of Jurkat cells. (A) Time courses for the directional migration of parental jurkat cells and **(B)** NuLight™ Red cells towards SDF-1α (50ng/ml) are comparable. In both studies, cells were plates at 5,000 cells per well in a 96-well IncuCyte™ ClearView™ plate and images were captured at 1 or 2 hour intervals in an IncuCyte® ZOOM system with 10x objective.

Optimization protocol

Unpacking and Storage Instructions

1. Check all containers and vials for leakage or breakage.
2. Remove the frozen vial from the dry ice packaging and transfer immediately to liquid nitrogen storage until ready for use.

Thawing and Culturing Cells

Jurkat NuLight™ Red Cells are supplied as non-adherent, immortalized T lymphocyte cells and should be maintained at a cell density of between 0.1-1 x 10⁶ cells/ml, do not allow the cell density to exceed 3 x 10⁶ cells/ml. Cell doubling time is approximately 20 hours. We recommend thawing one vial of Jurkat NuLight™ Red cells into one T25 flask with 6 mL of media.

1. Remove the vial of Jurkat NuLight™ Red cells from liquid nitrogen storage.
2. Thaw the vial by gentle agitation in a 37°C water bath. Be careful not to submerge entire vial to avoid contamination.

This process should take no more than 2 minutes. Remove the vial when only a tiny ice crystal remains.

3. Wipe the vial with 70% ethanol.
4. Transfer the vial contents to a centrifuge tube containing 9 mL of growth media and centrifuge at 125 x g for 5 min.
5. Remove the supernatant taking care not to disturb the cell pellet.
6. Resuspend the cells in 6 mL of growth media and add to a T25 flask

Note: To maintain expression of the NuLight™ Red label, we recommend culturing the cells in growth media containing 0.5 µg/ml puromycin. Puromycin can be removed for experiment/assay set-up.

7. Incubate the cells at 37°C incubator, 5% CO₂.
8. Cultures can be maintained by the addition or replacement of fresh growth medium.

Licenses and Warranty

Essen BioScience warrants that this product performs as described on the product label and in all literature associated with the sale of said product when used in accordance with the described protocol. This limited warranty is the sole warranty. No other warranties exist that extend beyond this warranty, either expressed or implied. Essen BioScience disclaims any implied warranty of merchantability or warranty of fitness for a particular purpose. Essen BioScience disclaims any responsibility for injury or damage and shall not be liable for any proximate, incidental or consequential damages in connection with this product.

If it is proven to the satisfaction of Essen BioScience that this product fails to meet performance specifications, Essen BioScience's sole obligation, at the option of Essen BioScience, is to replace the product or provide the purchaser with credit at or below the original purchase price. This limited warranty does not extend to anyone other than the purchaser. Notice of suboptimal performance

must be made to Essen BioScience within 30 days of receipt of the product.

This Essen BioScience product contain proprietary nucleic acid(s) coding for proprietary fluorescent protein(s) being, including its derivatives or modifications, the subject of pending patent applications and/or patents owned by Evrogen JSC (hereinafter "Evrogen Fluorescent Proteins"). The purchase of Essen BioScience products incorporating these fluorescent proteins conveys to the buyer the non-transferable right to use Evrogen Fluorescent Proteins only for research conducted by the buyer. No rights are conveyed to modify or clone the gene encoding fluorescent protein contained in this product or to use Evrogen Fluorescent Proteins for commercial purposes. The right to use Evrogen Fluorescent Proteins specifically excludes the right to validate or screen compounds for commercial purposes. For information on commercial licensing, contact Evrogen Licensing Department, email: license@evrogen.com.